




**LIGHTWEIGHT
AND FLEXIBLE**

Patslide Patient Transfer

The smooth upper surface provides low resistance, allowing for easy transfer. When the patients weight is applied, the non-skid inserts on the under surface grip firmly onto the operating table, trolley or bed. It's easy to position, and lightweight, weighing less than five kilos. The PATSLIDE even flexes to facilitate movement of the patient. Often siderails do not retract far enough causing transfer difficulty. The PATSLIDE provides a smooth bridge, spanning the gaps that occur between beds, trolleys and tables.

Features:

- Facilitates safe and easy patient transfers from bed to trolley etc
- Reduces the risk of back injury for the carer
- Handholes provide maximum grip for the carer
- Slip resistant pads increase stability during transfers
- Radiolucent, lightweight, impact resistant and anti-static
- Not to be used as a lifting device



PRODUCT CODE	PRODUCT NAME	DIMENSIONS
TSMS-LB0310	Patslide - All Purpose	1525mm (L) x 635mm (W)
TSMS-LB0320	Patslide - Theatre	1525mm (L) x 560mm (W)

JDHG100820-v1



To place a purchase order, contact us today:
sales@jdhealthcare.com.au | 1300 791 404

